

# Diagnostic Accessories

## SafeSEAL™ Antimicrobial Soft Diaphragm



Doctors Research Group, Inc.

SafeSEAL Antimicrobial Soft Diaphragm is impregnated with an Antimicrobial Compound, which preserves the integrity of the diaphragm, and thus makes it suppress the growth of bacteria and mold. As an improvement over traditional hard diaphragms, SafeSEAL is made from the highly-sensitive, yet durable C-Flex non-latex synthetic elastomer. The diaphragm conforms to the diameter of your stethoscope chestpieces to create an airtight seal, similar to that of a drum. SafeSEAL's Antimicrobial Soft Diaphragm is intended for use to minimize surface contamination and facilitate cleaning and disinfection of the stethoscope diaphragm. Standard precautions should be followed, including cleaning of diaphragms between use. SafeSEAL's unique construction is also resistant to ripping or hardening over time.



- |             |  |                                |
|-------------|--|--------------------------------|
| SS25-A-VAR  | SafeSEAL Adult Variety Pack  | (Contains 6 of each color + 1) |
| SS25-P-VAR  | SafeSeal Pediatric Variety Pack  | (Contains 6 of each color + 1) |
| SS50-4CL-4C | SafeSEAL individually foil package in a 4 column, 4 color dispenser box for weekly changing, Color-coded calendar on side of box | (50 of each standard colors)   |

Request a sample  
1-888-883-7331

Standard Colors: Green, Blue, Orange, & Red

The Spectrum® Dual Head and Spectrum® Nurse stethoscope features chrome-plated brass binaural, lightweight anodized aluminum chestpiece, 22" vinyl Y-tubing, spare diaphragm and a pair of mushroom ear tips. Overall length is 30".

### Spectrum® Dual Head

- 4601 Blue
- 4602 Black
- 4607 Burgundy
- 4616 Teal



Other colors available upon request

### Spectrum® Nurse

- 4281 Blue
- 4282 Black
- 4287 Burgundy
- 4286 Teal



Other colors available upon request

## Precision™ Series Aneroid Sphygmomanometer

Reusable with Calibrated Blue Nylon Cuff - Complete Kit with case

The Precision™ Series is ideal for the cost conscious healthcare provider who is looking for quality and affordability. This series features a durable cuff with hook and loop closure. The gauge is backed by a ten-year calibration warranty. Includes a standard air release valve and zippered carrying case. Latex-Free models include a five year inflation system warranty.

- 6011 Adult
- 6013 Infant
- 6015 Child
- 6016 Large Adult
- 6017 Thigh

**LATEX FREE**

- 6111 Adult
- 6113 Infant
- 6115 Child
- 6116 Large Adult
- 6117 Thigh



## Pulse Time Teaching Stethoscope - 4 Parties

Do you hear what I hear? No longer do medical instructors have to describe the faint sounds to students. With the Teaching Stethoscope, instructor and students listen to the same source of sound that offers multi-functional teaching applications and the necessary flexibility in all clinical and educational settings.

13-757PT Navy Blue



## Mobile Aneroid Sphygmomanometer

The Mobile Aneroid Sphygmomanometer incorporates all the features of the Desk and Wall Aneroid Sphygmomanometer on a quiet, sturdy, easy-rolling stand. Perfectly suited for vicinities in hospitals, clinics & nursing homes where repeated blood pressure measurements are conducted. 100% latex free

13-830MS



## Desk & Wall Aneroid Sphygmomanometer

The super-duty Desk & Wall Aneroid Sphygmomanometer system features calibrated manometer for clinical grade, certified accuracy. The safety, durability and convenience features make it ideal for high traffic areas in the emergency room, clinic or ambulance.

13-840DW

