

Breathing Bags and Accessories

Test Lung

Our test lung is dipped from blue Neoprene with a cemented external band. It is supplied with a standard 22mm neck insert and is available with or without foam inserts.

- Item #675AP Test Lung with out foam insert
- Item #676AP Test Lung with foam insert



RUSCH

Disposable Breathing Bags Green Neoprene



- ▶ Disposable
 - ▶ 22mm I.D. connection
 - ▶ Latex-Free Green Neoprene
- (A) 6201 1 liter Breathing Bag
 - (B) 6202 2 liter Breathing Bag
 - (C) 6203 3 liter Breathing Bag

Reusable Breathing Bags Blue Neoprene



RUSCH

- ▶ Reusable
- ▶ 22mm I.D. connection
- ▶ Latex-Free Blue Neoprene

- 6-6150 500cc Breathing Bag with Rolled Neck
- 6-6155 500cc Breathing Bag
- 6-6100 1 liter Breathing Bag
- 6-6102 2 liter Breathing Bag
- 6-6103 3 liter Breathing Bag
- 6-6104 4 liter Breathing Bag
- 6-6105 5 liter Breathing Bag

RUSCH

Reusable Breathing Bags

Our Fluted Breathing Bags are manufactured of black conductive neoprene and are produced by the dipping process. Bags are supplied with a tail and a loop for hanging. These Fluted Bags, often called pecan-shaped or football-shaped by our customers, are longer and narrower than a traditional oval bag and will provide an easier grip for smaller hands.

- 6-2170 Pecan Shaped Neoprene Breathing Bag 2 liter X 22 mm
- 6-2175 Pecan Shaped Neoprene Breathing Bag 1 liter X 22 mm
- 6-2193 Pecan Shaped Neoprene Breathing Bag without Neck Insert 500cc
- 6-2194 Pecan Shaped Neoprene Breathing Bag 500cc X 22mm
- 6-2200 Pecan Shaped Neoprene Breathing Bag 1500cc X 22mm

