

# CO2 Easy, Resuscitator Set, and Manujet III


## CO2 Easy

Featuring the largest visual indicator available, the CO2 Easy is ideally mounted up front where clinicians and first responders can effectively view critical carbon dioxide levels without looking away from the patient. The non-toxic, pH sensitive indicator uses familiar purple-to-yellow color conversion when CO2 is exhaled.

- ▶ Designed for fast, easy airway assessment and verifies proper endotracheal tube placement.
- ▶ Helps prevent unwarranted intubation and confirms cardiac output during CPR.
- ▶ Expansive 2" diameter Visual Indicator makes rapid breath-to-breath response easy to see.
- ▶ Anti-fogging agent reduces susceptibility to moisture and humidity, ensuring continuous lens clarity.
- ▶ Complies with Standard of Care Guidelines from the American Heart Association, American College of Emergency Physicians, American Association for Respiratory Care and the American Society of Anesthesiologists.

8-2134 CO2 Easy Adult  
8-2135 CO2 Easy Pediatric

	<u>Adult Co2 Easy</u>	<u>Pediatric CO2 Easy</u>
Internal Volume	30 cc	3cc
Resistance to Flow	3.0cm H2O +/- 1cm at 60 L/min	2.5cm H2O +/- 0.5 cm at 10 L/min
Weight	Less than 20g	Less than 5g
Connector Ports	Patient End: 22mm OD/15mm ID Proximal End: 15mm OD	Patient End: 18mm OD/15mm ID Proximal End: 15mm OD/5mm ID
Patient Use	33 lbs. or above	2.2-33 lbs or less



## Resuscitator Set with Case

- ▶ Complete Set Includes: Resuscitation Bag, Mask, Reservoir Bag, Reservoir Valve, 2mm Oxygen Tubing
- ▶ Latex Free
- ▶ PVC Set is Disposable
- ▶ Silicone Set is Reusable - Autoclaveable up to 134°C

Items can be purchased individually



**VBM**

## Disposable PVC Resuscitator Set

- 8-8519 Adult: 1800ml Resuscitation Bag, Mask Size 5, Reservoir Bag 2600ml, Reservoir Valve, and 2m O<sub>2</sub> Tubing
- 8-8529 Child: 500ml Resuscitation Bag, Mask Size 3, Reservoir Bag 2600ml, Reservoir Valve, and 2m O<sub>2</sub> Tubing
- 8-8539 Infant: 250ml Resuscitation Bag, Mask Size 1, Reservoir Bag 600ml, Reservoir Valve, and 2m O<sub>2</sub> Tubing

## Reusable Silicone Resuscitator Set

- 8-8019 Adult: 2000ml Resuscitation Bag, Mask Size 5, Reservoir Bag 2600ml, Reservoir Valve, and 2m O<sub>2</sub> Tubing
- 8-8029 Child: 500ml Resuscitation Bag, Mask Size 3, Reservoir Bag 2600ml, Reservoir Valve, and 2m O<sub>2</sub> Tubing
- 8-8039 Infant: 250ml Resuscitation Bag, Mask Size 1, Reservoir Bag 600ml, Reservoir Valve, and 2m O<sub>2</sub> Tubing

## Jet Ventilation Manujet III

- ▶ The Manujet III is supplied in a case with 4m (13ft.) pressure hose, Jet-Ventilation Catheters acc. to Ravussin (13G, 14G, 16G), 100 cm (3.3ft.) connecting tube, Endojet Adapter (for ET tube) and Bronchoscope Adapter
- ▶ Latex Free
- ▶ For Difficult Airway Management
- ▶ Ideal for EMS and emergency carts

Item #8-1003 Manujet Set



**RUSCH**