

## B & F Oxygen Pressure Regulators For Quality That Costs Less!




B&F regulators offer an array of options, such as an easy to read flow rate window and a dial-in click-style flow controller ... making operation convenient for the variety of users. B&F oxygen regulators are available with or without a hose barb outlet, with or without a DISS outlet, and with a one-year manufacturer's warranty.

- ▶ Uni-body construction designed to strengthen the regulator and prevent breakage.
- ▶ Nickel-plated brass body for safe use in oxygen service.
- ▶ Two body designs: Yoke inlet for small cylinders, nut and nipple inlet for large cylinders.
- ▶ Compact, 1-1/2" diameter low-profile contents gauge with gauge protector.
- ▶ Serrated dial-in flow control knob for ease of flow rate selection.
- ▶ Large easy to read flow rate numbers.
- ▶ Hose barb or DISS outlet provide secure connection of oxygen supply line to patient.



**CGA 870 Yoke with Hose Barb**

10-21020	0-6 LPM
10-21014	0-8 LPM
10-21017	0-15 LPM
10-21030	0-25 LPM



**CGA 540 Nut**

10-23020	0-6 LPM, Hose Barb
10-23002	0-8 LPM, DISS outlet
10-23003	0-15 LPM, DISS outlet

**CGA 540 Nut with Hose Barb and Check Valves**

10-23004	0-15 LPM, 1 valve
10-23006	0-15 LPM, 2 valves
10-23025G	0-25 LPM, 2 valves

**CGA 870 Yoke with DISS**

10-22020	0-6 LPM
10-23014	0-8 LPM
10-23017	0-15 LPM
10-21031	0-25 LPM

**CGA 870 Yoke with Hose Barb and Check Valves**

10-21012	0-15 LPM, 2 valves
10-24017	0-15 LPM, 1 valve
10-21025	0-25 LPM, 2 valves
10-21032	0-25 LPM, 1 valve

### Neonatal and Pediatric Liquid Oxygen Flowmeter

Model 320, 3 LPM oxygen flowmeter calibrated for 20 to 25 psi liquid oxygen, is intended for use with neonatal or pediatric liquid oxygen systems that deliver 20 to 25 psi oxygen. Supplied with 1/8" NPT female back inlet and a DISS male oxygen bottom outlet.

Item # 10-0210 Neonatal or Pediatric Liquid Oxygen Flowmeter 3 LPM

### Cylinder Wrench for D & E Cylinders

- 7480 Metal wrench
- 7480-C Metal wrench with security chain
- 7479 Plastic wrench pk/10
- 7479-C Plastic wrench with security chain

