

ET Tube Stylets & Exchangers

METTS (Muallem ET Tube Stylet)

for difficult intubation
for single use, sterile

core made of metal, malleable (maintains curvature)
flexible pre-formed tip, graduation marks



Item #	Size	Length	for ET Tube
7-3208	8CH	40 cm	>/-3.5 mm
7-3212	12 CH	65 cm	>/-5.0 mm
7-3214	14 CH	65 cm	>/-6.0 mm

Introducer

for difficult intubation
for single use, sterile

Lumen-design, rigid tip, graduation marks



Item #	Size	Length	for ET Tube
7-3414	14 CH	65 cm	>/-5.0 mm

Stylets

Short, malleable intubation stylet
autoclavable up to 134 degrees C

Soft, atraumatic tip



Item #	Size	Length	for ET Tube
Without adjustable 15 mm ID connector			
7-3635	small	34 cm	2.5 - 3.5 mm
7-3645	medium	38 cm	4.0 - 4.5 mm
7-3665	large	45 cm	5.0 - 6.5 mm
7-3611	large	45 cm	7.0 - 11.0 mm

Item #	Size	Length	for ET Tube
With adjustable 15 mm ID connector			
7-3735	small	34 cm	2.5 - 3.5 mm
7-3745	medium	38 cm	4.0 - 4.5 mm
7-3765	large	45 cm	5.0 - 6.5 mm
7-3711	large	45 cm	7.0 - 11.0 mm

METTI (Muallem ET Tube Introducer)

for difficult intubation and tube exchange
for single use, sterile

core made of plastic, semi-rigid (flexible)
flexible pre-formed tip, graduation marks



Item #	Size	Length	for ET Tube
7-3312	12 CH	80 cm	>/-5.0 mm
7-3314	14 CH	80 cm	>/-6.0 mm

Tube Exchanger

for extubation and exchange of Endotracheal Tubes
for single use, sterile

Lumen-design with distal side holes allows
oxygenation via included 15 mm or Luer Lock
Connector, graduation marks




Item #	Size	Length	for ET Tube
7-3511	11 CH	80 cm	>/-4.0 mm
7-3514	14 CH	80 cm	>/-5.0 mm
7-3519	19 CH	80 cm	>/-7.0 mm

Tube Exchanger

for extubation and exchange of Endotracheal Tubes
for single use, non-sterile

Smooth, Rounded Tip, High Density Polyethylene,
graduated marks

A simple, rapid technique for changing endotracheal
tubes. Can be done in less than two minutes!



Stylettes & Endotracheal Tube Changers

Disposable Malleable Aluminum Stylet



Item#	Length cm	Length inches	Size O.D.
3650	40cm	16"	8Fr
3655	30cm	11.5"	8Fr

New Lower Price

Flexi-Slip Stylet

RUSCH



- ▶ Sterile
- ▶ Plastic coated malleable metal core
- ▶ Color coded
- ▶ Coated for easy withdrawal
- ▶ Single Use

Item#	Length	Size O.D.
3-2501	30cm	6 Fr
3-2503	44cm	10 Fr
3-2505	45cm	12 Fr
3-2507	45cm	14 Fr

New Lower Price

Slick® Disposable Endotracheal Stylet

RUSCH



- ▶ Sterile
- ▶ Nylon-coated wire
- ▶ Malleable
- ▶ Sureloc® Stop helps maintain desired length for patient safety
- ▶ Handle forming capability for better grip
- ▶ Single Use

Item#	Length	Size O.D.
3-3106	32cm	6 Fr
3-3108	42cm	8 Fr
3-3110	45cm	10 Fr

Sheridan T.T.X. Exchanger

- ▶ Allows quick and efficient endotracheal tube exchange or replacement
- ▶ Flexible for easy insertion
- ▶ Frosted surface minimizes drag
- ▶ Includes depth marks
- ▶ Hollow to allow for spontaneous breathing
- ▶ Sterile



3-2402	2.0 - 4.0mm
3-2404	4.0 - 6.0mm
3-2406	6.0 - 8.5mm
3-2408	8.5 - 10.0mm

HUDSON RCI

Sheridan Jet Ventilation/ Tracheal Tube Exchanger

- ▶ Jet ventilation capability with the JETTXX provides additional time to assess spontaneous breathing
- ▶ Luer lock adapter on JETTXX can be easily attached to conventional jet ventilation devices
- ▶ Can be used for oxygen delivery
- ▶ For use with tube size 6.5 - 10.0mm ID
- ▶ Depth marks provide a guide for placement of endotracheal tube, without use of a laryngoscope
- ▶ Hollow to allow for spontaneous breathing
- ▶ Sterile



3-2420 6.5 - 10.0mm