


The Timeter® Classic & Soft-Touch Flowmeters



Neonates to adults, the entire series of time-tested Timeter® Classic Flowmeters demonstrates our unsurpassed knowledge of the needs of care providers and patients alike. You can see it in the accuracy. You can see it in the craftsmanship. You can see it in the performance, year after year. Because we understand that special care with extreme accuracy is critical to fragile, little bodies, Classic Neonatal Flowmeters are specially designed. When you're looking for a superior, long-lasting flowmeter that's easy to maintain and offers extreme accuracy, you'll find it in our Timeter® Classic series.



Model Max 100 psi	Oxygen Flowmeter 200 cc	Oxygen Flowmeter 1 LPM	Oxygen Flowmeter 3 LPM	Oxygen Flowmeter 16 LPM	Oxygen Flowmeter 75 LPM	Air Flowmeter 16 LPM	Air Flowmeter 75 LPM
Back Fitting	Extremely Low Flow	Very Low Flow 1/10-1 LPM	Low Flow 1/8 - 3 LPM	General Purpose 1/4-16 LPM	High Flow IMV & C-PAP 5-75 LPM	General Purpose ½-16 LPM	High Flow IMV & C-PAP 5-75 LPM
3/4" Fitting 1/8" x 1/8"	10-T01118	10-T01211	10-T0060	10-T0010	10-T0030	10-T0020	10-T0040
Chemetron Male	10-T03118	10-T03221	10-T0360	10-T0310	10-T0330	10-T1620	10-1640
Oxequip Male	n/a	n/a	10-T0160	n/a	10-T0130	10-T4020	n/a
Ohmeda Male	10-T04118	10-T04211	10-T0460	10-T0410	10-T0430	10-T1920	10-1940
Puritan Bennett	10-T02118	10-T02211	10-T0260	10-T0210	10-T0230	n/a	n/a
DISS Female HT	10-T12118	10-T12211	10-T1260	10-T1210	10-T1230	n/a	10-T2040
DISS Female	10-T06118	10-T06211	10-T0660	10-T0610	10-T0630	n/a	n/a
DISS Male	n/a	10-T10211	10-T1060	10-T1010	10-T1030	n/a	10-1540

The Soft-Touch Timeter® Flowmeter is durable, long-lasting and easy to maintain. Quality plus Value. Flowmeter hood is made of crack-resistant polycarbonate. Easier readability on the 0-5 LPM range. Steel float eliminates sticking caused by static electricity. Nickel-plated solid brass body provides durability and long life. Larger, easier to grip Soft-Touch knob. Color coded for air or oxygen.



Back Fitting	Oxygen Flowmeter 15 LPM	Oxygen Flowmeter 15 LPM Power Take-off	Oxygen Flowmeter 8 LPM	Oxygen Flowmeter 8 LPM Power Take-off	Air Flowmeter 15 LPM	Air Flowmeter 15 LPM Power Take-off
3/4" Fitting 1/8" x 1/8"	10-15103	10-15103P	10-08103	10-08103P	10-15113	10-15113P
Chemetron Male	10-15203	10-15203P	10-08203	10-08203P	10-15213	10-15213P
Oxequip Male	10-15303	10-15303P	10-08303	10-08303P	10-15313	10-15313P
Ohmeda Male	10-15403	10-15403P	10-08403	10-08403P	10-15413	10-15413P
Puritan Bennett	10-15603	10-15603P	10-08603	10-08603P	10-15613	10-15613P
DISS Female HT	10-15803	10-15803P	10-08803	10-08803P	10-15813	10-15813P
DISS Female	10-15903	10-15903P	10-08903	10-08903P	10-15913	10-15913P
DISS Male	10-15003	10-15003P	10-08003	10-08003P	10-15013	10-15013P