

Flowmeters and Neb Block



Flowmeter



- ▶ Chrome body flowmeter,
- ▶ Black glass floats indicator ball
- ▶ 5 year warranty
- ▶ Low Flow, Pediatric, and Compact Models Available

Oxygen Flowmeter 1-15 LPM

- 7501-PM 1/8" NPT female inlet
- 7503-PM DISS Handtight Nut
- 7507-PM Ohmeda male inlet
- 7509-PM Chemetron male inlet
- 7511-PM Puritan male inlet

Air Flowmeter 1-15 LPM

- 7501-APM 1/8" NPT female inlet
- 7503-APM DISS Handtight Nut
- 7507-APM Ohmeda male inlet
- 7509-APM Chemetron male inlet
- 7511-APM Puritan male inlet



Dial Flowmeter

- ▶ For headwall outlets or ceiling mount drops
- ▶ Two different flow ranges: 1/5-15 or 1-25 LPM
- ▶ Power take-off and Air Dial Flowmeters available



- 7503-DF15 1/2 -15 LPM DISS Handtight Nut
- 7507-DF15 1/2 -15 LPM Ohmeda male inlet
- 7509-DF15 1/2 -15 LPM Chemetron male inlet

- 7503-DF25 1-25 LPM DISS Handtight Nut
- 7507-DF25 1-25 LPM Ohmeda male inlet
- 7509-DF25 1-25 LPM Chemetron male inlet

For connection on Top, add "T" to end of Part Number.

DISS, Oxequip, Puritan-Bennett & Schrader Available.

Oxygen Select Flowmeter



- ▶ 0-15 LPM, Lightweight, Metal DISS male connection
- ▶ Easy to use. Simply push the button on the bottom of the unit and rotate the desired port
- ▶ Lifetime warranty and the needle valves carry a five year warranty
- ▶ Eliminates tubing nipple adaptors (it has 2 of its own)

- 7501-SF 1/8" NPT female inlet
- 7503-SF DISS Handtight Nut
- 7507-SF Ohmeda male inlet
- 7509-SF Chemetron male inlet
- 7511-SF Puritan male inlet
- 7513-SF Oxequip male inlet



Flow Selector

- ▶ Give nebulizer treatments without disconnecting the humidifier
- ▶ Reduces changeover time to seconds
- ▶ No more loose tubing laying around
- ▶ Eliminates tubing nipple adaptors (it has 2 of its own)
- ▶ Fits all flowmeters that have DISS male outlet fittings
- ▶ Lightweight

7519 Flow Selector

The Flow Selector is a selectable single flow device. Flow is only delivered to the selected connection indicated by the green "on" arrow.

Neb Block

It Increases Productivity ...

Neb treatments, resuscitation bags or equipment power are all available with this unique product. All equipment can be kept connected and ready for use when needed.

It's Versatile ...

Need power for equipment? The Neb Block is equipped with a 50 psi power take off to allow use of any equipment that requires 50 psi power.

It Saves Time ...

You simply need to push the switch from "OFF" to "ON" to deliver a neb treatment or utilize a resuscitation bag.

It's Convenient ...

The Neb Block is available with either chrome or compact flowmeters. Connections for all wall outlet styles are available.

It's Preset ...

The outlet barb is available in 6 pre-set flows for easy delivery of treatment.



DISS Handtight Inlet with Chrome Flowmeter

- 7503-NB1 6-7 LPM
- 7503-NB2 7-8 LPM
- 7503-NB3 8-9 LPM
- 7503-NB4 10-11 LPM
- 7503-NB5 14-16 LPM
- 7503-NB6 24-26 LPM

Ohmeda Inlet with Chrome Flowmeter

- 7507-NB1 6-7 LPM
- 7507-NB2 7-8 LPM
- 7507-NB3 8-9 LPM
- 7507-NB4 10-11 LPM
- 7507-NB5 14-16 LPM
- 7507-NB6 24-26 LPM

Chemetron Inlet with Chrome Flowmeter

- 7509-NB1 6-7 LPM
- 7509-NB2 7-8 LPM
- 7509-NB3 8-9 LPM
- 7509-NB4 10-11 LPM
- 7509-NB5 14-16 LPM
- 7509-NB6 24-26 LPM