

Gas Sampling Lines, Nasal Mask, Cannulas, & Oxygen Tubing

Gas Sampling Lines & Filter



- ▶ Manufactured to the highest standards of quality from medical grade material to fit all standard luer connections.
- ▶ We offer PVC lines for end tidalCO2 monitoring and coextruded lines for anesthetic agent monitoring, which are to be used when measuring anesthetic gas concentrations.

7735	PVC/PE Coextrusion, .050" ID X 10', Male X Male
7736	PVC/PE Coextrusion, .050" ID X 10', Male X Female
7737	PVC, Male X Male, .060" ID X 10', Male X Male
7739	Microphobic Filters, Clear, Male / Female Luers

BiFlo™ & PediFlo™ Nasal Mask

Westmed



- ▶ The BiFlo nasal mask is a simple and comfortable device designed as an alternative for oxygen therapy patients who cannot or will not tolerate a cannula or for those patients who have recurring necrosis.
- ▶ PaO2 and FIO2 values are comparable to those receiving oxygen with a conventional nasal cannula.
- ▶ The BiFlo delivers a gentle flow of oxygen. Unlike a nasal cannula, oxygen flow is not jetted into the nasal cavity.
- ▶ The PediFlo animated design reduces the apprehension generally encountered by younger patients.

8-B106	PediFlo Bunny nasal mask with 7' Kink Resistant Tubing
8-B101	BiFlo Adult nasal mask with 7' Kink Resistant Tubing
8-B125	BiFlo Adult nasal mask with 25' Kink Resistant Tubing

LATEX FREE

Comfort Plus™ Soft Tipped Cannula

Westmed

- ▶ With the patient's ease of fit and maximum comfort in mind, the Comfort Plus line of cannulas are designed lightweight and clear to ensure a comfortable fit for long-term use.
- ▶ The Over-the-Ear, "Bolo-Style" adjustable slide moves smoothly, yet provides a secure and comfortable fit while allowing maximum freedom of movement.
- ▶ The curved nasal prongs enhance and promote alignment of the cannula to fit naturally on the patient's face, eliminating the irritating flap required with conventional style cannulas.
- ▶ Comes with standard "six channel" safety tubing to ensure uninterrupted flow of oxygen to the patient.



7-1637	Adult Micro w/ 7' Kink Resistant Tubing
7-1673	Infant w/ 7' Kink Resistant Tubing and Threaded Nut
7-1683	Pediatric w/ 7' Kink Resistant Tubing and Threaded Nut
7-1687	Pediatric w/ 7' Kink Resistant Tubing
7-1690	Adult Harness only w/ 2" Connector
7-1697	Adult w/ 7' Kink Resistant Tubing
7-1698	Adult w/ 14' Kink Resistant Tubing & long straight conn
7-1699	Adult w/ 14' Kink Resistant Tubing and Threaded Nut
Conventional Cannulas	
7-1610	Adult Harness w/ 2" connector
7-1617	Adult with 7' Kink Resistant Tubing
7-1614	Adult w/ 14' Kink Resistant Tubing

LATEX FREE

Oxygen Tubing

Westmed



- ▶ Kink resistant six channel inner lumen design
 - ▶ Available in two durometer options:
75 shore (clear) and 85 shore (green)
- | | |
|--------|--|
| 7-0007 | 7' Kink Resistant Tubing |
| 7-0008 | 7' Kink Resistant Tubing w/ Threaded Nut |
| 7-0014 | 14' Kink Resistant Tubing |
| 7-0027 | 25' Kink Resistant Tubing (Green) |

LATEX FREE