


Medical Gas Wall Outlets, Swivel Hose Adaptor, and Manifold

Medical Gas Wall Outlets

Available in Four Keying Styles for Most Medical Gases



X

Quick Connect Type:
 0 = non (rough in only)
 1 = Ohmeda
 2 = Chemetron
 3 = Puritan Bennett
 4 = DISS

Y

Gas Service:
 1 = Vacuum
 2 = Medical Air
 3 = WAGD
 4 = Carbon Dioxide
 7 = Nitrogen
 8 = Nitrous Oxide
 9 = Oxygen

Z

Configuration:
 0 = none (latch only)
 1 = 90 degree tubing
 2 = 180 degree tubing
 3 = Hose barb

How to order: Prefix is MWS

Wall outlet with DISS (X) + Nitrous Oxide (Y) + 90 degree tubing (Z) = Item #MWS-481

Please call for a quotation, 1-888-883-7331

90° Swivel Hose Adaptor



This 90° Swivel hose adaptor is revolutionizing hose fittings in the medical field and beyond. The innovative design provides increased mobility while maintaining a secure connection; and you can save money on replacements. The swivel adaptor helps prevent dangerous kinks and chaffing that can result in leakage.

90° Swivel Hose Adaptor with 1/4" Conductive 6 Foot Hose

Oxygen

HA-S51072H1H5 Ohmeda Male X Ohmeda Female
 HA-S51072C1C5 Chemetron Male X Chemetron Female

Air

HA-S31072H2H6 Ohmeda Male X Ohmeda Female
 HA-S31072C4C6 Chemetron Male X Chemetron Female

Suction

HA-S41072H3H7 Ohmeda Male X Ohmeda Female
 HA-S41072C3C7 Chemetron Male X Chemetron


Swivel Only

8-FPT24 1/4", 90° Swivel Adaptor X Hose Barb
 8-FPT25 5/16", 90° Swivel Adaptor X Hose Barb

Available for all medical gases, adaptors, quick connects, and lengths. Please call 1-888-833-7331.

Extension Manifold - DISS M/Check Valve

- ▶ Economical expansion of single wall outlet
- ▶ Permanently mounts to wall
- ▶ 2 foot conductive hose with 1/8" NPT Female inlet is standard



10-7230 Oxygen 2 outlet
 10-7231 Oxygen 3 outlet
 10-7232 Oxygen 4 outlet

10-7233 Air 2 outlet
 10-7234 Air 3 outlet
 10-7235 Air 4 outlet