

Gas Sampling Lines, Nasal Mask, Cannulas, & Oxygen Tubing

Gas Sampling Lines & Filter




- ▶ Manufactured to the highest standards of quality from medical grade material to fit all standard luer connections.
- ▶ We offer PVC lines for end tidal CO₂ monitoring and coextruded lines for anesthetic agent monitoring, which are to be used when measuring anesthetic gas concentrations.

7735-25 PVC/PE Coextrusion, .050" ID X 10', Male X Male
7735 PVC/PE Coextrusion, .050" ID X 10', Male X Male
7739 Microphobic Filters, Clear, Male / Female Luers

BiFlo™ & PediFlo™ Nasal Mask

Westmed



- ▶ The BiFlo nasal mask is a simple and comfortable device designed as an alternative for oxygen therapy patients who cannot or will not tolerate a cannula or for those patients who have recurring necrosis.
- ▶ PaO₂ and FIO₂ values are comparable to those receiving oxygen with a conventional nasal cannula.
- ▶ The BiFlo delivers a gentle flow of oxygen. Unlike a nasal cannula, oxygen flow is not jetted into the nasal cavity.
- ▶ The PediFlo animated design reduces the apprehension generally encountered by younger patients.

LATEX FREE

8-B106 PediFlo Bunny nasal mask with 7' Kink Resistant Tubing
8-B101 BiFlo Adult nasal mask with 7' Kink Resistant Tubing
8-B125 BiFlo Adult nasal mask with 25' Kink Resistant Tubing

Comfort Plus™ Soft Tipped Cannula

Westmed

- ▶ With the patient's ease of fit and maximum comfort in mind, the Comfort Plus line of cannulas are designed lightweight and clear to ensure a comfortable fit for long-term use.
- ▶ The Over-the-Ear, "Bolo-Style" adjustable slide moves smoothly, yet provides a secure and comfortable fit while allowing maximum freedom of movement.
- ▶ The curved nasal prongs enhance and promote alignment of the cannula to fit naturally on the patient's face, eliminating the irritating flap required with conventional style cannulas. Standard "six channel" safety tubing to ensure uninterrupted flow.

7-1652 Premature with 7' Kink Resistant Tubing
7-1677 Infant with 7' Kink Resistant Tubing
7-1687 Pediatric with 7' Kink Resistant Tubing
7-1637 Adult Micro with 7' Kink Resistant Tubing
7-1697 Adult with 7' Kink Resistant Tubing
7-1696 Adult with 14' Kink Resistant Tubing
7-1699 Adult with 14' Kink Resistant Tubing and Threaded Nut

Specialty Cannulas

7-2503 Adult CO₂/O₂ cannula with male luer
7-2504 Adult CO₂/O₂ cannula with female luer
7-2507 Adult Demand Cannula



LATEX FREE

Oxygen Tubing

Westmed

- ▶ Kink resistant six channel inner lumen design
- ▶ Available in two durometer options:
75 shore (clear) and 85 shore (green)

7-0007 7' Kink Resistant Tubing
7-0008 7' Kink Resistant Tubing w/ Threaded Nut
7-0014 14' Kink Resistant Tubing
7-0027 25' Kink Resistant Tubing (Green)
7-0040 40' Kink Resistant Tubing
7-0050 50' Kink Resistant Tubing

LATEX FREE