


# Nasopharyngeal Airways and The Easy Tube

## NasoSafe Latex Free Nasal Airway

NasoSafe is a new and innovative Nasopharyngeal Airway. The swivel safety grip eliminates the risk of the airway from traveling down the nasal passage while providing the maximum comfort for the patient.



- ▶ The NasoSafe swivel grip protects the Nasopharyngeal airway from loss in the nasal passage.
- ▶ The swivel grip allows for positioning of the airway in the most comfortable position.
- ▶ The swivel safety grip is color coded to match suction catheter size color coding.
- ▶ The swivel safety grip can be positioned so that it does not block the free nostril.
- ▶ The rounded edge of the tip combined with angle cut design helps to reduce the risk of injury during intubation.
- ▶ Latex Free




1238-60	NasoSafe Nasal Airway 6.0 ID Green
1238-70	NasoSafe Nasal Airway 7.0 ID Orange
1238-80	NasoSafe Nasal Airway 8.0 ID Red
1238-90	NasoSafe Nasal Airway 9.0 ID Purple

### RUSCH

### Robertazzi Latex Free Nasal Airway

- ▶ Manufactured from soft neoprene
- ▶ Smooth beveled tip for gentle insertion
- ▶ Resilient
- ▶ Trumpet end for secure placement
- ▶ Sterile, single patient use
- ▶ 9 sizes



Item#	Size
1234-20	20 Fr
1234-22	22 Fr
1234-24	24 Fr
1234-26	26 Fr
1234-28	28 Fr
1234-30	30 Fr
1234-32	32 Fr
1234-34	34 Fr
1234-36	36 Fr

Nasal Airway Kits

Item #1234-06 Sizes 22 - 32 Fr  
 Item #1234-09 Sizes 20 - 36 Fr

### RUSCH

### The EasyTube

Compare to Nellcor Combitube

Designed for difficult tracheal or esophageal intubation. The EasyTube creates a viable airway whether it is placed in the trachea or in the esophagus. When placed in the esophagus, the EasyTube allows the use of fiber-optic devices, or passage of a suction catheter or tube exchanger. The risk of tracheal trauma is reduced to a minimum due to the single lumen at the distal tip.



- ▶ Single lumen at distal tip reduces potential for trauma
- ▶ High volume low pressure PVC cuff at distal end allows for minimal wall pressure
- ▶ Soft Latex Free cuff to seal the pharynx
- ▶ Special placement rings clearly indicate a visible depth marking that can be used for either tracheal or esophageal intubation
- ▶ Large volume syringe, pre-filled with air for the pharyngeal cuff, and small volume syringe, pre-filled with air for the distal cuff, making intubation quick and easy

1841-41	41 Fr	Large
1828-28	28 Fr	Small

The complete EasyTube set includes tube, syringes pre-filled with air, and an extra length suction catheter. Packaged in a convenient roll-up kit.