

Oxygen Analyzers

Oxygen Analyzer/Monitor



The MAXO₂® is known for its legendary 2-year warranty on both the analyzer and sensor. Provides simple operation and exceptional stability. Long battery life and overall low cost of ownership.

- Item #8-OM25ME MAXO₂® Oxygen Monitor
(Includes External O₂ Sensor Hi/Low Alarms)
- Item #8-OM25AE MAXO₂® Oxygen Analyzer
(Includes External O₂ Sensor)
- Item #8-OM25A MAXO₂® Oxygen Analyzer
(Includes Internal O₂ Sensor)

Replacement Sensors:
 Item #8-MAX250E Sensor for 8-OM25ME & 8-OM25AE Monitors (24 month warranty)
 Item #8-MAX250 Sensor for 8-OM25A Monitors (24 month warranty)

The Handi® Oxygen Analyzer - 2 Year Warranty

The most unique oxygen analyzer of its class and size. It can be used anywhere from the NICU to homecare applications.

- ▶ Handi will last between 2 - 4 years
- ▶ Pocket Size
- ▶ Simple, reliable, affordable
- ▶ 3 digit LCD display
- ▶ 2 year warranty



Item #8-207P01 The Handi® with Lanyard

Oxygen Analyzer

The Pacifitech Oxygen Analyzer is designed to provide continued or spot-checking of oxygen delivery devices in respiratory care and oxygen therapy applications. The Analyzer is available in an internal or external sensor mounting configuration. All analyzers are supplied with sensors.

Instrument Features: Oversize LCD Digital Display, Low Battery Warning, Electro-Magnetic Compatibility (EMC) with EN 6-601-1-2, Belt Clip (internal sensor version), Plastic wedge for pole or wall mounting (external sensor version) Compact Size, One Year Warranty
 Sensor Features: Galvanic O₂ partial pressure, Maintenance free, External or internal sensor options, Temperature compensated, One Year Warranty

Oxygen Analyzer Specifications:

Warm Up Time: None required
 Display Range: 0-100% Oxygen
 Display Resolution: 0.1% Increments
 Accuracy: +/- 2% Full Scale
 Response Time: <10 seconds to 90% final value
 Sensor Life: 1 year
 Power Supply: 9-volt Alkaline Battery
 Battery Life: Approximately 3000-4500 hours
 Weight: 7.1 oz. (200 G)
 Temp Range: 0°C to 40°C (32°F to 104°F)



PT1000E External Mounted Sensor
 PT1000I Internal Mounted Sensor

Replacement Sensors:
 PT08A Sensor for External Sensor
 PT08I Sensor of Internal Sensor

MAXO2+ Oxygen Analyzer and Replacement Oxygen Sensors

MAXO2+ Oxygen Analyzer



The revolutionary MAXO2+ oxygen analyzer is fully equipped with simple to use one-touch calibration. Both analyzer and sensor are protected with a two-year warranty and a unsurpassed 5000 hour battery-life. Patented design is user friendly and fits comfortably in your hand.

8-OA25A Oxygen Analyzer with Internal Sensor
 8-OA25AE Oxygen Analyzer with External Sensor

Replacement Sensor
 8-MAX250P For 8-OA25A Internal Sensor
 8-MAX250E For 8-OA25AE External Sensor



Replacement Sensors for Anesthesia and Respiratory Applications

New Lower Price



Anesthesia

Item#	Description	Warranty
8-001	Datex Ohmeda Excel Modulus 4700 Oxicap, 5250 RGM, 5120/5125 Monitor, Hudson 5556, Mercury 10-103-07, Teledyne R15, VTI 101272 & 100352	(12 month warranty)
8-011	Dual Cathode Draeger Narkomed Incubators, IC8000, Hudson 5557, Mercury 10-103-07, Teledyne R23MED, VTI 101372 & 100352	(12 month warranty)
8-008	Datex Ohmeda Excel Modulus 5250 RGM ext only, 5120/5125 Monitor, Mercury 10-103-13, Teledyne R32MED, VTI 101272	(18 month warranty)
8-002	Narkomed Series 6400, Mercury 10-103-12, Teledyne R31MED, & VTI 101372	(18 month warranty)
8-010	Datex Ohmeda 7900 Series Ventilator, Mercury 10-103-10	(12 month warranty)
8-22C	Criticare Systems 1100 Series Ventilator & Poet Models, Teledyne R22MED	(12 month warranty)



Respiratory

Item#	Description	Warranty
8-009	Hudson 5500, Mercury 10-103-03, Teledyne R13	(14 month warranty)
8-012	Hudson 5566, Mercury 10-103-06, Teledyne R24MED, Seimens 9004979-E347E	(14 month warranty)
8-013	Hudson 5567, Mercury 10-103-05, MiniOx I,II,III & MiniOX 3000, Mallinckrodt 7820 Monitor & 7200 Ventilator, Teledyne R17MED, & VTI 103100	(14 month warranty)
8-014	MiniOX I-A (Internal Sensor)	(14 month warranty)
8-016	Mallinckrodt Ventilators, Puritan Bennett	(12 month warranty)
8-017	Hudson 5558, Mercury 10-103-02, Teledyne T-7	(12 month warranty)
8-018	Hudson 5568, Mercury 10-103-01, Teledyne C1R	(12 month warranty)
8-019	Hudson 5569, Mercury 10-103-00, Teledyne T1	(12 month warranty)